

# Quazar DM 6000

Instruction Material for Quazar DM 6000  
Precision O.E.M. Laser Diode Module and Driver Unit

# Quazar

## Diode Module and Driver

### Replacement Parts

#### Laser Supplies

#### Laser Eyewear

Eyewear # D-213-4600 Each \$199.99

Laser Module Replacement 6000 mW

6 Watt Laser Module # D-213-6700 Each \$599.99

Carbon Dye 50ml Bottle

Carbon Dye # P-213-2500 Each \$29.99

Prices are subject to change without notification. To order on-line, go to  
<http://www.centre-biotechnique-avance.com>  
For technical assistance beyond what this manual provides, please e-mail  
[admin@centre-biotechnique-avance.com](mailto:admin@centre-biotechnique-avance.com)  
Please allow 24 hours for processing.

## Quick Setup Guide

Read this guidebook first to set up your  
equipment for use.



Keep this manual in a convenient place for quick and easy  
reference at all times.

The product names in this guidebook are trademarks or registered trade marks of each specific manufacturer. In the interest of providing superior equipment, Quazar Industries reserves the right to modify or amend equipment specifications without notice or obligation.

## Important Information

- Read all safety and operating instructions before connecting or using this equipment.
- To protect against electrical shock, do not use this equipment near water. Do not immerse unit, plug transformer while standing in water, or spray with any liquids. As with many electrical appliances, this unit carries a live charge even while unplugged. Do not dismantle this unit (there are hazardous voltages inside).
- Do not place this unit near an open flame or cooking/heating devices (e.g., stoves, heat registers, radiators, etc.).
- If an extension cord is used, it should be appropriately rated for voltage, power, and frequency as indicated on the back of the unit.
- The power cord should be routed so that it is not likely to be walked upon or pinched (especially near the wall outlet, extension receptacle, or where the cord exits the unit).
- To avoid serious damage to the power system and microelectronics, the transformer should be unplugged from the wall outlet when the unit will be unused for long periods of time.

## Trouble Shooting

Should you encounter technical problems with your Quazar DM 6000 Photo Epilator, refer to the following guide for potential problems and their solutions.

**—Unit is plugged into the wall, all accessories are correctly inserted into the unit but no laser output is being registered.**

- ++**Check all connections. Plug and unplug each one being sure all contacts are sound.**
- ++**Check all cords. Due to continual bending and fatigue, wires may fray or break resulting in full loss of power.**
- ++**Check Fuse:** The Quazar Epilator unit has a fast-acting fuse mounted on the top panel of your unit to protect the delicate laser diode from voltage spikes, on-line power surges and electrostatic discharge (ESD damage). Replace with 400-500 micro amp (ua) fast acting fuse only. Failure to comply with these specifications may result in serious damage to your laser and will void all warranties.

**—Unit hums or makes noises.**


- ++**Unit needs servicing.**
- Laser output is weak.**
- ++**Unit needs servicing.**

**—No output from the laser is registered after all trouble-shooting suggestions listed above have been checked.**

- ++**Unit needs servicing.**

**—Laser beam is not focused to a usable point (output is a weak line or oval).**

- ++**Your optics need to be adjusted. Normally, the lens is pre-set by the manufacturer to focus at roughly 2-3 inches from the aperture. If the user has altered this setting (turned the optical spanner adjustment) without correct instructions the beam may be unusable until corrected. Turn the lens adjuster with a straight-slot screw driver (see notches at right) to reset beam focus. Turning clockwise will bring the focus out (making the point farther away). Turning counterclockwise will bring the point closer.**



## Equipment Warranty

We warrant to the original purchaser the equipment manufactured by us to be free from defects in material and workmanship under normal use and service. Our obligation under this warranty shall be limited to the repair or exchange of any part or parts which may prove defective under normal use and service within 12 calendar months from the date of shipment and which our examination shall disclose to our satisfaction to be thus defective. When necessary, purchaser shall apply for a Return Materials Authorization and instructions on proper return procedures from their original sales associate. The laser diode (head) requires special operating precautions which, if defied, may void warranty.

## Warranty Extension Certification:


Customer Number \_\_\_\_\_ Authorization Number \_\_\_\_\_  
Warranty Extension ( ) years Warranty Type: A B C D

## The Treatment Procedure

Shine the red dot on the highlighted hair follicle. Adjust the distance of the laser head from the tissue to create a pin-point of focus. The optimal beam diameter for maximum intensity is roughly 1 mm.

The follicle will begin to flash as the photon energy reacts with the carbon dye. It is also normal to see some gaseous emissions (vapor and smoke). Continue to hold the laser in place until all visible reactivity ceases. This may take up to 15 seconds. Move on to the next follicle and repeat.

Some patients may find the laser treatment uncomfortable. In these cases the use of a topical anesthetic such as lidocaine (a synthetic amide, C14H22N2O, used chiefly in the form of its hydrochloride as a local anesthetic and antiarrhythmic agent) will reduce discomfort. A mild cryogenic (nitrogen-based) topical spray will also minimize any pain.



**Treatment Around or Near the Eyes:** Great care must be exercised when working near the eyes. The laser emission is powerful enough to actually penetrate the eyelid and permanently damage the eye. Having the patient close their eyes is not satisfactory protection. The use of a dark-colored damp wash cloth which is folded over four times will deflect the harmful radiation; however, only laser-protective eyewear is recommended.

**Treatment Around or Near Mucus Membranes:** Laser radiation will severely damage the tissues inside the nose and ear canal. Treatment should be avoided in these areas altogether.

**Treatment Around or Near the Genitals:** Laser hair removal is safe for application to the pubic regions including the reproductive organs or both sexes. Care must be taken into consideration in these areas due to the increased level of neural sensitivity. The patient may find the process uncomfortable without a topical.

**Treatment Around or Near the Areola (nipple):** Laser hair removal is safe and effective on hair growth which occurs from the areola of both sexes. Again, care must be taken into consideration in these areas due to the increased level of neural sensitivity.

**Post-Treatment:** The skin surrounding the treatment area may experience short-term erythema (reddening) which will subside within 12 hours. Should the treatment area show signs of excess scabbing you may wish to reduce the overall treatment time or intensity. The application of a post-treatment cooling and healing gel (such as Aloe) is encouraged to speed healing and reduce sensitivity. Instruct the patient to refrain from applying cosmetics or sunbathing for at least 24 hours.

## Treatment Efficacy

Permanent hair removal is a gradual process which takes 90 days or more for complete destruction of the follicle tissues. Each hair must go through its entire growth cycle for it to be effectively treated. Only during the **early anagen** phase is it vulnerable to destruction. The following **chart** will give you an accurate example of what the reduction in growth activity should look like from 30, 60, and 90 days of treatments.



**Pointers:** Human skin has roughly 1,000 follicles per square inch. Only a small percentage of them are active (anagen) at any given time. For effective treatment, the hair needs to be treated during the anagen phase. It normally takes about 90 days for each follicle to go through the full growth cycle. Should you see hair growth activity in an area which has been completely cleared, it is most likely from a follicle which was dormant at the time of the original treatment. Simply apply treatment to these hairs as they appear. All hair growth activity will normally terminate within 3 months.



Patient #	Hair counts		6 week clearance percent
	pre	post	
1 Females, 3 males	274	9	97%
	327	19	94%
	3	2	99%

The first laser hair removal treatment was administered, quite by accident, by a technician in the late 1970's. He was repairing a high output 692 nm YAG laser when he unwittingly left his arm in the beam's path for a few seconds. Although his skin was completely unharmed, all the hair in the area had completely burned off. In the years to follow, the area of his arm which was exposed to the laser remained completely **bold**.

## Laser 'Electrolysis'

What actually occurred was a "heat exchange" reaction with the pigment (called melanin) deep inside the follicle tissue. Photon energy from the laser had penetrated into the translucent dermis virtually unobstructed. When it reached the pigment naturally contained in the follicle, radiation quickly heated those cells (*melanocytes*) to well over 100 degrees Celsius. Human tissue cannot survive at this temperature, henceforth, the complete destruction of the follicle was virtually instantaneous. The scientific term for this process is called **thermolysis**.



Laser light with the wavelength of 600 to 900 nm (nanometer) passes through human tissue with very little loss of intensity. The Quazar DMD500 Epilator produces a precise wavelength of 808 nm, which is proven to provide the greatest tissue penetration while limiting the loss of energy to natural pigments found in the skin.

Quazar also utilizes a constant-wave output laser (CW), which makes it one of the most effective and sophisticated systems in the world. Other systems operate strictly on a 'pulsed-type' output. The problem with pulsed lasers is that follicle tissues are only heated for a very small increment of time (usually less than 1/1,000,000 of a second). Recent clinical studies show that it takes a full **two seconds** or more for complete carbonization, desiccation and coagulation at 100 degrees Celsius. With Quazar, the beam may be held in position until full destruction is achieved. By carefully regulating the output at 808 nm, the laser will not harm the skin, even after several minutes of constant exposure.

The best candidate for laser hair removal has fair skin with dark terminal hairs. Skin typing for exposure to ultraviolet light can be categorized by the Fitzpatrick classification, developed by Dr. Thomas Fitzpatrick of Harvard Medical School.

**Skin Type I:** Never tans, always burns (extremely fair skin, blonde hair, blue/green eyes)

**Skin Type II:** Occasionally tans, usually burns (fair skin, sandy to brown hair, green/brown eyes)

**Skin Type III:** Often tans, sometimes burns (medium skin, brown hair, brown eyes)

**Skin Type IV:** Always tan, never burns (olive skin, brown/black hair, dark brown/black eyes)

**Skin Type V:** Never burns (dark brown skin, black hair, black eyes)

**Skin Type VI:** (black skin, black hair, black eyes)

Types I through 4 are outstanding candidates. Type 5 will have excellent results as well but care must be taken to assure that the laser will not burn the skin. This is achieved by using a cryogenic spray or air-flow accessory. Type 6 should not undergo laser hair removal unless used in conjunction with skin bleaching due to the high risk of burning and hypo/hyper pigmentation issues.

- References**
1. Bjerring P, Carranos M, Egekvist H, Christiansen K, Trullius A. Hair reduction using a new intense pulsed light irradiator and a normal mode ruby laser. *J Clean Laser* 2000; 1: 3-7.
  2. Anderson RR. Treatment of pseudofolliculitis with a pulsed infrared laser. *Arch Dermatol* 2000; 136:1343-6.
  3. Emswile S, Li CC, Newman N. Laser hair removal with alexandrite versus diode laser using four treatment sessions: 1-year results. *Dermatol Surg* 2001; 27: 925-9.
  4. Gargul N, Aslan G, Akce T, Erdogan B. Comparison of alexandrite laser and electrolysis for hair removal. *Dermatol Surg* 2000; 26:37-41.
  5. Benayahu N, Luid A, Galimberti M, Ferretti G. Long-term ablation with long-pulsed red/infrared: YAG laser. *Dermatol Surg* 1999; 25:175-8.
  6. Bercini JK, Mihov M. Long-term evaluation of the long-pulsed alexandrite laser for the removal of thin hair at shortened treatment intervals. *Dermatol Surg* 2000; 26:633-7.


## Warnings and Advisories

**Quazar DM 6000 Epilator** produces laser radiation which can be harmful to the eyes. Always wear protective eyewear while operating this equipment. Laser radiation has the capability to burn the skin if the technician does not closely observe the patient's reaction to the procedure.

**Advisory:** Laser electrolysis results in full destruction of the hair follicle and is **irreversible**. Always plan ahead before undertaking detail work such as eyebrow shaping or hairline contouring. Patch test a small area (no larger than 1X1 inch square) before full application. Allow 24 hours to determine the patient's reaction.

This manual provides a tutorial overview of galvanic electrolysis. For more detailed information, please refer to "**Modern Electrolysis, Volume 2, Phototherapy**" by Palamed Press. If you are unsure of the proper use of this device, do not use. The **Quazar DM 6000 Epilator** is intended for use by qualified individuals or professionals only.

## Control Locations/Feature Descriptions



- A. **Key Lockout:** This feature is required by law on all high-power laser devices. The first step in the correct sequence to power-up your laser is to turn this lock clockwise using the special key included in your kit. The light (center LED lamp) will be green when the power system is off and red when the power system is on.
- B. **System Status LED:** This light will indicate system status (neutral or live). Green is neutral and red is live.
- C. **Laser Jack:** A four-prong coaxial power jack for the diode laser module.
- D. **Power Cord:** Rated for 60Hz, 120-240 V, 10 Amp with corresponding plug-style for country of destination.
- E. **Laser Optics:** This feature allows the technician to control the size, shape and characteristic of the laser beam. Using the optics spanner wrench, turn clockwise to make the beam focus at a shorter distance, and counterclockwise to increase the focal length. It is set at 1.5 inches at the factory.
- F. **Eyewear:** This is an essential part of the treatment process. Direct or reflective laser radiation can seriously injure the eye. Both the technician and the patient must use the protective eyewear while the laser is enabled. Eyewear is intended for **accidental** exposure only. Never stare directly into a laser beam.
- G. **Carbon Dye:** This is an 'atomized' form of molecular carbon which easily penetrates deeply into the follicle shaft. The dye adds pigment which gives a receptor for the photon/heat exchange reaction. The carbon atoms will capture the laser energy and convert it into heat for the rapid and efficient cauterization of tissue for the permanent destruction of the hair follicle organ.
- H. **High-Precision Tweezers:** Apparatus for the extraction of follicle prior to carbon dye application.
- I. **Trigger Switch:** This red button should be activated by way of the operator's thumb. It will initiate laser pulses from the laser module.



**Warning:** Laser radiation is emitted from this aperture. Unprotected eye exposure may cause serious injury resulting in loss of vision or blindness. Always use laser eyewear.

## Pre-Treatment

Before applying treatment, remove all hair from the area by tweezing or waxing. Laser hair removal is most effective when applied to an empty follicle shaft. Human hair simply does not have enough pigment to allow for sufficient heat exchange to cauterize, desiccate and necrotize the cells which produce hair. To compensate for this lack of 'quantitative' and 'qualitative' photon targets, it will be necessary to place a high-density carbon dye inside the follicle prior to treatment.



will block the dye. The hairs may be left in the skin if desired, but results will be improved if they are extracted.

Isolate the hairs to be permanently destroyed. It will be necessary to remove them by swiftly plucking (or waxing) in the direction of growth. Pulling slowly generally leaves most of the follicle tissue inside the pore which results will be improved if they

## Photo-Reactive Dye Application

Using a cotton-tipped applicator, completely cover the treatment area with the special dye included in your kit. Massage the dye into the follicle pore with a firm downward circular motion. Repeat 2-3 times to saturate the follicle pore. Use an **ethyl alcohol** based wipe (isopropyl alcohol will not dissolve the dye) to lightly clean the excess from the surface of the skin. At this point you will have all desired follicles **visibly highlighted** with a dark spot (as seen above) and are ready to power up your laser for treatment. Carbon dye must be also used if the hair is left in the skin for laser 'shaving'.



**Pointer:** If your patient objects to having a depilatory process before the treatment, you may continue without the carbon dye. This alternate procedure is the equivalent of 'laser shaving' (for which long-term permanency is marginal). For best results, the use of a photo-reactive dye is highly recommended.

## Dermal Coolant Application

Place a thin layer of laser dermal coolant spray on the treatment area prior to laser application. This will protect the surface of the skin from burning as well as improve the translucency of the skin (rate at which light can pass). Failure to use the dermal coolant prep may result in unnecessary discomfort for the patient during treatment and increase the likelihood of a surface burn. Should the liquid become dry, it will be necessary to re-apply frequently. The use of a humidifier in dry climates will substantially prolong the duration for which the dermal prep will retain its cooling properties.

## Laser Startup Procedure

1. Before plugging your laser head into the power system, be sure the key switch is off (LED will show green). **Failure to adhere to the start-up procedure may cause irreversible damage to your laser head (which would not be covered by warranty).**
2. Plug the laser module accessory into the power system, tighten the spanner nut completely.
3. Place the laser head in a secure position **POINTED AWAY** from you and the patient.
4. Power up your module by using your key lockout switch (turn to the 'enabled' setting).
5. An audible-tone will accompany the activation of the laser. This is an additional safety feature to warn the operator that the module is producing radiation.



Your laser is now **live** and extreme care must be exercised while handling to avoid accidental eye-exposure. The output of the laser is invisible (infrared) at 808 nm. Only 1% of the photon emission is made visible for the purpose of placement and/or focusing.

You must turn off your laser in **reverse order** of the start-up sequence. Failure to do so may cause irreversible damage to your laser head and is not covered by warranty. Always use caution when turning on and/or operating this system. The laser is very powerful and can cause **serious eye injury** if used incorrectly.